

STOP B&W INODORE CON INDICATORE

ODOURLESS B&W STOP WITH pH INDICATOR

BWINDEXSTOP1

BWINDEXSTOP500

DESCRIZIONE DESCRIPTION

Bagno di arresto inodore con Indicatore di pH per il trattamento Bianco e Nero da usare per fermare l'azione dello sviluppo sia nel processo Carta che Pellicola.

Il prodotto è disponibile nel formato da 500 mL e da 1 litro di concentrato da diluire 1 + 19.

Odorless stop bath with pH Indicator dye for Black and White paper and film processing. Use to stop development. Available in a 500 mL bottle or a 1 litre bottle of concentrate.

ISTRUZIONI DI SICUREZZA SAFETY AND HANDLING INSTRUCTIONS

Indossare guanti, occhiali ed indumenti protettivi. Tenere fuori dalla portata dei bambini.

Fare riferimento alla scheda di sicurezza del prodotto, MSDS o contattare la Bellini Foto al +39 075 985 174.

For detailed safety and handling instructions refer to MSDS or contact BELLINI FOTO S.r.l. by e-mail:

info@bellinifoto.it or by telephone +39 075 985 174. Handle with care and follow safety and handling instructions.

Wear gloves and use safety goggles; keep out of reach of children.

SPECIFICHE RACCOMANDATE SPECIFICATIONS

Diluizioni: 1 + 19 Tempo di trattamento: 10 - 30 secondi a 20 - 25 °C

Dilution: 1 + 19 Processing time: 10-30 seconds at 25 °C

ISTRUZIONI MISCELAZIONI MIXING INSTRUCTIONS

	Acqua Water	Concentrato Concentrate
1 L	950 mL	50 mL

